

Gojira Fine Chemicals, LLC

Quality Products & Services Meeting Customer Needs
5386 Majestic Parkway, Suite 7, Bedford Heights, OH 44146
www.qojirafc.com • sales@qojirafc.com

Specification

Lysing Matrix E, 2 ml Tubes

Catalog Number: LM1007

CAS Number: 6914-010; 6914-050; 6914-100; 6914-500

Description: Lysing Matrix E, distinguished by purple cap is optimized for DNA and

RNA extractions from complex matrices such as: soil, geological sediments. It also consists of environmental samples such as:

wastewater and filtrates, forensic samples, and infected tissues. It is the

most cited in literature lysing matrix for processing of soil and environmental samples. It consists of optimized mixture of 0.1mm diameter borosilicate beads, 1.4 mm diameter zirconium silicate beads and 4mm diameter glass beads, all packed in an impact resistant tube.

Date: June 4, 2021

Storage: 15° - 30°C

Test	Specification
Physical Description	A Mixture of White Fine Beads, Off-White Ceramic
	Spheres and 4 mm Glass Beads in a 2 ml micro-centrifuge
	tube.
Identity Test	Passes
DNase	Not Detected
RNase	Not Detected

Approved: John Grubb

John Grubb, PhD Gojira Fine Chemicals

Quality Manager

Tel: 440-252-5397, Fax: 888-211-5523

docsupport@gojirafc.com